



Point of Care Informatics

# POCcelerator® System Device Interface List

Total number of interfaces: 230

Last updated: April 1, 2023

[siemens-healthineers.com/informatics](https://siemens-healthineers.com/informatics)

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
4S Information Systems	DAWN	Supplier software	→			
Abaxis	Piccolo Xpress <b>UPDATE</b>	Clinical chemistry	→			
Abbott Diabetes Care	FreeStyle Precision Pro (FSPP) <b>UPDATE</b>	Glucose, Ketone	↔	•	•	•
	FreeStyle Precision Pro (FSPP) Wireless mode <b>UPDATE</b>	Glucose, Ketone	↔	•	•	•
	PCx Precision	Glucose	→			
	Precision Xceed Pro (PXP)	Glucose, Ketone	↔	•	•	•
Abbott Point of Care	Afinion AS100 <b>UPDATE</b>	HbA1c	↔	•		•
	Afinion 2 <b>UPDATE</b>	HbA1c	↔	•		•
	Alere i (ASTM)	Multitest	→			
	Alere i (POCT1-A) <b>UPDATE</b>	Multitest	↔	•	•	•
	DIGIVAL <b>UPDATE</b>	Multitest	→			
	CELL-DYN Emerald	Hematology	→			
	Cholestech LDx	Multitest	→			
	ID NOW <b>UPDATE</b>	Multitest	↔	•	•	•
	Pima	HIV	→			
	i-STAT 1 (via i-STAT-DE-DMI software) <b>UPDATE</b>	Multitest	↔	•	•	•
i-STAT Alinity <b>UPDATE</b>	Multitest	↔	•	•	•	
Acon Laboratories	Mission U120	Urinalysis	→			
Aidian (former Orion Diagnostica)	QuikRead 101	CRP	→			
	QuikRead go (ASTM)	CRP	→			
	QuikRead go (POCT1-A2) <b>UPDATE</b>	CRP	↔	•	•	•
A. MENARINI Diagnostics	Zenit RA	Immunoassay	→			
Analyticon	CombiScan 100	Urinalysis	→			
	Urilyzer 100	Urinalysis	↔	•		
	Urilyzer 100 Pro	Urinalysis	↔		•	
	Urilyzer 100 Pro V2 / DocUReader 2 Pro <b>UPDATE</b>	Urinalysis	↔	•	•	•
Arkray	Adams HA-8160	HbA1c	→			

→ Unidirectional connectivity    ↔ Bidirectional connectivity

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Beckman Coulter	Coulter Ac-T diff	Hematology	→			
	Coulter Ac-T diff 2	Hematology	→			
	Coulter Ac-T 8	Hematology	→			
	Coulter Ac-T 10	Hematology	→			
	DxH 500	Hematology	→			
	DxH 520	Hematology	→			
Becton Dickinson	BD Veritor Plus (via BD Cloud Connectivity Solution)	Multitest	→			
Biomérieux	mini VIDAS	Multitest	→			
Bio-Rad	Bio-Rad Unity Real Time (URT)	Supplier software	→			
Biosynex	BSX Reader	Multitest	→			
biosurfit	Spinit	Multitest	→			
Boditech	AFIAS-1	Immunoassay	→			
	AFIAS-6	Immunoassay	→			
Bohle Medical	Medonic M-series M16/M20	Hematology	→			
	Medonic M-series M16C/M20C	Hematology	→			
	Medonic M-series M16M/M20M	Hematology	→			
	Medonic M-series M16S/M20S	Hematology	→			
	Medonic M-series M32B	Hematology	→			
	Medonic M-series M32C	Hematology	→			
	Medonic M-series M32M	Hematology	→			
	Medonic M-series M32S	Hematology	→			
	Swelab Alfa Basic / Plus Basic	Hematology	→			
	Swelab Alfa Standard / Plus Standard	Hematology	→			
	Swelab Alfa Cap / Plus Cap	Hematology	→			
	Swelab Alfa Sampler / Plus Sampler	Hematology	→			
Caretaker Medical	Caretaker4	Vital Signs Monitor	→			
Cepheid	GeneXpert PC and Omni	Multitest	→			
	GeneXpert Infinity	Multitest	→			
	GeneXpert Xpress	Multitest	↔	•		•
Chirus	DXpress Reader	Multitest	→			
CoaguSense	Coag-Sense PT2	Coagulation	↔	•	•	
concile	Omega 100 (Ω100)	Immunoassay	→			
Diagnosics For The Real World (DRW)	Samba II	Multitest	↔		•	
Diatron	Diatron Aquila	Hematology	→			
DnaNudge	NudgeBox	Immunoassay	↔	•		
Dräger	Jaundice Meter JM-105	Bilirubin	→			
Drucker Diagnostics	QBC Star	Hematology	→			
EKF Diagnostics	Quo-Lab (via CIB) <b>UPDATE</b>	HbA1c	↔	•	•	•
	Quo-Test (via CIB) <b>UPDATE</b>	HbA1c	↔	•	•	•
Enigma Diagnostics	Enigma MiniLab	Multitest	→			
ET Healthcare	ET Pylon	Multitest	→			

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Eurolyser Diagnostica	Eurolyser CUBE	Multitest	→			
	Eurolyser CUBE-S	Multitest	→			
	Eurolyser smart	Multitest	→			
	Eurolyser solo	Multitest	→			
Fuji	Fuji Dri-Chem 4000	Clinical chemistry	→			
	Fuji Dri-Chem NX10N	Clinical chemistry	→			
	Fuji Dri-Chem NX500	Clinical chemistry	→			
	Fuji Dri-Chem NX600	Clinical chemistry	→			
	Fuji Dri-Chem NX700	Clinical chemistry	→			
Genedrive	Genedrive System <b>UPDATE</b>	Multitest	→			
GenePOC Diagnostics	Revogene <b>UPDATE</b>	Multitest	↔	•		•
Haemonetics	ClotPro <b>UPDATE</b>	Hemostasis	↔		•	
	TEG 6s (via TEGManager HaemoCommunicator)	Hemostasis	↔		•	
Helena Laboratories	Cascade Abrazo <b>UPDATE</b>	Coagulation	↔			•
	Cascade POC	Coagulation	↔	•		
HemoCue	Glucose 201 DM	Glucose	↔	•	•	•
	Glucose 201 DM RT	Glucose	↔	•	•	•
	Hb 201 DM	Hemoglobin	↔	•	•	•
	HbA1c 501	Hemoglobin	→			
	WBC Diff	Hematology	→			
HemoSonics	Quantra	Hemostasis	↔	•	•	
HLS Therapeutics	CSAN Pronto	Hematology	→			
Hologic	PeriLynx	Immunoassay	→			
	Rapid fFN	Immunoassay	→			
HORIBA	ABX Micros 60	Hematology	→			
	ABX Micros ES 60	Hematology	→			
	ABX Microsemi CRP LC-667 <b>NEW</b>	Hematology	→			
	ABX Microsemi CRP LC-767 <b>UPDATE</b>	Hematology	→			
	Yumizen H500 <b>UPDATE</b>	Hematology	→			
	Yumizen H550 <b>NEW</b>	Hematology	→			
Lite-On Technology Corp.	skyla HB1 (POCT1-A2) <b>UPDATE</b>	Clinical chemistry	↔	•		
	skyla HB1 (ASTM)	Clinical chemistry	→			
	skyla Hi	HbA1c	↔	•		
LumiraDx	LumiraDx (Instrument and Host) <b>UPDATE</b>	Multitest	↔	•	•	•
Macherey-Nagel	Uryxxon Relax (Connect) <b>UPDATE</b>	Urinalysis	↔		•	
Magnasense Technologies	Magnia Reader	Multitest	→			
Medtronic	ACT Plus	Coagulation	→			
	HMS Plus	Coagulation	→			
MeMed	MeMed Key	Multitest	→			
Mitsubishi	PATHFATH	Immunoassay	→			
nal von minden (BioSciTec and Mediwiss)	Rapid Slide Scanner (RSS) III	Drug-Screen	→			

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Nihon Kohden	Celltac Alpha (MEK 6500/6510)	Hematology	→			
	Celltac Chemi (CHM-4100)	Hematology	→			
Norma Diagnostika	Norma Icon 3	Hematology	→			
	Norma Icon 5	Hematology	→			
Nova Biomedical	NovaNet*	Supplier software	↔	•	•	
	Stat Profile pHox Ultra (via NovaNet)	Blood gas	→			
	Stat Profile Prime (via NovaNet)	Blood gas	→			
	Stat Profile Prime Plus (via NovaNet)	Blood gas	↔	•	•	
	StatSensor Creatinine (via NovaNet)	Creatinine	↔	•	•	
	StatStrip Glucose (via NovaNet)	Glucose	↔	•	•	
	StatStrip Glucose/Ketone (via NovaNet)	Glucose, Ketone	↔	•	•	
	StatStrip Hb/Hct (via NovaNet)	Hemoglobin	↔	•	•	
	StatStrip HHL (via NovaNet)	Hemoglobin/ Hematocrit/Lactate	↔	•	•	
	StatStrip Lactate (via NovaNet)	Lactate	↔	•		
	Stat Profile CCX	Blood gas	↔	•	•	
	Stat Profile pHox / pHox plus	Blood gas	↔		•	
	Stat Profile pHox Ultra	Blood gas	→			
	Stat Profile Prime	Blood gas	→			
	StatSensor Creatinine <b>UPDATE</b>	Creatinine	↔	•	•	•
	StatStrip Glucose <b>UPDATE</b>	Glucose	↔	•	•	•
	StatStrip Glucose/Ketone <b>UPDATE</b>	Glucose, Ketone	↔	•	•	•
	StatStrip Hb/Hct <b>UPDATE</b>	Hemoglobin	↔	•	•	•
	StatStrip Lactate <b>UPDATE</b>	Lactate	↔	•	•	•
	OPTI Medical	OPTI CCA-TS2 <b>UPDATE</b>	Blood gas	↔	•	•
Philips	BiliChek	Bilirubin	→			
PixCell Medical	HemoScreen <b>UPDATE</b>	Hematology	↔	•	•	•
Quidel	Sofia	Immunoassay	→			
	Sofia 2	Immunoassay	→			
	Triage Meter Plus	Cardiac markers	↔	•		
	Triage Meter Pro	Cardiac markers	↔	•		
Radiometer	ABL7xx series	Blood gas	↔		•	
	ABL8xx series	Blood gas	↔		•	
	ABL80 FLEX series	Blood gas	↔		•	
	ABL835 (POCT1-A) <b>NEW</b>	Blood gas	→			
	ABL90 FLEX series	Blood gas	↔		•	
	ABL90 (POCT1-A) <b>NEW</b>	Blood gas	→			
	AQT90 FLEX	Multitest	↔		•	
	AQUIRE	Supplier software	↔		•	
	NPT7	Blood gas	→			
RADIANCE	Supplier software	↔		•		
Response Biomedical	RAMP 200 <b>UPDATE</b>	Multitest	↔	•		•

\*NovaNet is needed for the StatStrip-Family connections with a firmware version greater than 4.17.8.47.  
Other devices can be connected directly or via NovaNet

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Roche Diagnostics	ACCU-CHEK Inform II <b>UPDATE</b>	Glucose	↔	•	•	•
	CARDIAC reader	Cardiac markers	↔			•
	CoaguChek Pro II <b>UPDATE</b>	Coagulation	↔	•	•	•
	CoaguChek XS Plus <b>UPDATE</b>	Coagulation	↔	•	•	•
	CoaguChek XS Pro <b>UPDATE</b>	Coagulation	↔	•	•	•
	cobas b 101	Hemoglobin	↔	•	•	•
	cobas b 123 <b>UPDATE</b>	Blood gas	↔	•	•	
	cobas c 111	HbA1c	→			
	cobas e 411	Immunoassay	→			
	cobas h 232 <b>UPDATE</b>	Cardiac markers	↔	•	•	•
	cobas 6000	Clinical chemistry	→			
	cobas IT 1000 <b>UPDATE</b>	Supplier software	↔	•	•	
	cobas Liat <b>UPDATE</b>	Multitest	↔	•		
	cobas u 411 <b>NEW</b>	Urinalysis	→			
	Reflotron Plus	Clinical chemistry	→			
	Reflotron Sprint	Clinical chemistry	→			
Urisys 1100	Urinalysis	↔	•			
Samsung	LABGEO <sup>HC10</sup>	Hematology	↔	•		
	LABGEO <sup>IB10</sup>	Immunoassay	↔	•		
	LABGEO <sup>PT10</sup>	Clinical chemistry	↔	•		
SD Biosensor	Standard F200	Multitest	→			
	Standard F2400	Multitest	→			

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
<b>Siemens Healthcare Diagnostics</b>	ADVIA® 360 Hematology System	Hematology	→			
	ADVIA Centaur® CP Immunoassay System	Immunoassay	→			
	ADVIA Centaur XP Immunoassay System	Immunoassay	→			
	Atellica® DCA Analyzer (Engine and IDM) <b>UPDATE</b>	Hemoglobin	↔	•	•	•
	Atellica 1500 Automated Urinalysis System	Urinalysis	↔			•
	Atellica UAS 800 Automated Urinalysis System	Urinalysis	↔			•
	Atellica VTLi Immunoassay System <b>UPDATE</b>	Immunoassay	↔	•		
	CLINITEK Advantus® Urine Chemistry Analyzer	Urinalysis	→			
	CLINITEK Status® Analyzer	Urinalysis	→			
	CLINITEK Status+ Urine Chemistry Analyzer <b>UPDATE</b>	Urinalysis	↔	•		•
	DCA Vantage® Analyzer	HbA1c	↔	•		•
	epoc® Blood Analysis System (Reader and Host) <b>UPDATE</b>	Multitest	↔	•	•	•
	epoc NXS Blood Analysis System (Reader and Host) <b>UPDATE</b>	Multitest	↔	•	•	•
	RAPIDComm® Data Management System	Supplier software	↔	•	•	
	RAPIDLab® 1200 series Blood Gas Systems <b>UPDATE</b>	Blood gas	↔	•	•	•
	RAPIDLab 348EX Blood Gas System <b>UPDATE</b>	Blood gas	→			
	RAPIDPoint® 500 Blood Gas System <b>UPDATE</b>	Blood gas	↔	•	•	
	RAPIDPoint 500e Blood Gas System <b>UPDATE</b>	Blood gas	↔	•	•	•
	Stratus® CS Acute Care™ Diagnostic System	Cardiac markers	→			
	Stratus CS 200 Acute Care™ Diagnostic System	Cardiac markers	↔	•		
<b>Sight Diagnostics</b>	OLO	Hematology	→			
	OLO (POCT1-A) <b>UPDATE</b>	Hematology	↔	•		•
<b>Sphere Medical</b>	Proxima	Blood gas	→			
<b>Spindiag</b>	Rhonda Player <b>UPDATE</b>	Multitest	↔	•		
<b>Stat-Diagnostica</b>	DiagCORE <b>UPDATE</b>	Multitest	→			
<b>Sysmex</b>	Alifax Roller	Hematology	→			
	CA-560	Coagulation	→			
	CA-600	Coagulation	→			
	pocH-100i	Hematology	→			
	UF-500i (via SUIT Manager)	Urinalysis	→			
	UF-1000i (via SUIT Manager)	Urinalysis	→			
	UX-2000 (via SUIT Manager)	Urinalysis	→			
	XE Series (via SUIT Manager)	Hematology	→			
	XN Series (via SUIT Manager)	Hematology	→			
	XN-L Series (via SUIT Manager)	Hematology	→			
	XS Series (via SUIT Manager)	Hematology	→			
	XT Series (via SUIT Manager)	Hematology	→			
	XP-300	Hematology	→			

Manufacturer	Device	Application	Connection	Download Operator	Positive patient ID	Configuration
Thermo Fisher Scientific	B-R-A-H-M-S Direct Reader <b>UPDATE</b>	Immunoassay	→			
Universal Biosensors	Xprecia Stride	Coagulation	↔	•		•
Veda.Lab	Easy Reader+ <b>UPDATE</b>	Multitest	↔	•		
Voden	VChemistry S	Clinical chemistry	→			
Welch Allyn	Connex Spot Monitor <b>UPDATE</b>	Vital Signs Monitor	→			
Werfen	Avoximeter 1000E	Blood gas	→			
	Avoximeter 4000	Blood gas	→			
	BIO-Flash	Immunoassay	→			
	GEM OPL <b>NEW</b>	Blood gas	→			
	GEM Premier 3000	Blood gas	↔		•	
	GEM Premier 3500	Blood gas	↔		•	
	GEM Premier 4000	Blood gas	↔	•	•	
	GEM Premier 5000	Blood gas	↔	•	•	
	GEMweb Plus <b>UPDATE</b>	Supplier software	↔	•	•	
	HEMOCHRON Response	Coagulation	→			
	HEMOCHRON Signature Elite <b>UPDATE</b>	Coagulation	↔	•		•
	ProTime InRhythm	Coagulation	→			
	ROTEM delta (direct or via Rotem live)	Hemostasis	↔		•	
	ROTEM gamma (direct or via Rotem live)	Hemostasis	↔		•	
	ROTEM sigma (direct or via Rotem live)	Hemostasis	↔		•	

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. By constantly bringing breakthrough innovations to market, we enable healthcare professionals to deliver high-quality care, leading to the best possible outcome for patients.

Our portfolio, spanning from in-vitro and in-vivo diagnostics to image-guided therapy and innovative cancer care, is crucial for clinical decision-making and treatment pathways. With our strengths in patient twinning, precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the biggest challenges in healthcare. We will continue to build on these strengths to help fight the world's most threatening diseases, improving the quality of outcomes, and enabling access to care.

We are a team of 66,000 highly dedicated employees across more than 70 countries passionately pushing the boundaries of what's possible in healthcare to help improve people's lives around the world.

Acute Care, ADVIA, ADVIA Centaur, Atellica, CLINITEK Advantus, CLINITEK Status, DCA Vantage, epoc, POCcelerator, RAPIDComm, RAPIDLab, RAPIDPoint, Stratus, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local representative for availability.

---

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

**Published by**

Siemens Healthcare Diagnostics Inc.  
Point of Care Diagnostics  
2 Edgewater Drive  
Norwood, MA 02062-4637  
USA  
Phone: +1 781-269-3000