

XL

Reliable support system for rehabilitation









#### Highlights

#### High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

#### Multifunctional

The treadmill can be used in various ergometry fields, like cardiology, pulmonary function, rehabilitation, sports medicine and physiotherapy, as a stand alone unit or controlled by other equipment.

#### Easy weight unloading

The body weight can be very quick and convenient unloaded by using the remote control

#### Mobile

The system can easily be moved when it is not used

#### Freedom of Movement

The test subject can move freely without any obstruction of the system







#### Reliable support system for rehabilitation

The Lode Body Weight Support System (BWSS) XL is especially developed to support patients that are not able to walk independently. They need support in order to be able to learn how to walk again. The Lode BWSS is the perfect product to enable this. During the training session, the patient weight can be supported accurately up to 80 kg, allowing a broad range of patients to use the system. The compensation weight can easily be adjusted with a remote control. The weight is compensated through a unique use of spring combinations, which is a breakthrough and provides the most stable compensation system. The XL version of the system is compatible with large surface treadmills (width > 60cm).

The height of the product is suitable for normal rooms (from 2,5 meters). For higher rooms you may consider to order height adjustment parts to enable use by taller patients. This will increase the height with 18 cm. Partnumber of this option is 932824.

#### **Features**



#### User friendly remote control

The user friendly remote control offers both adjustment possibilities through the Versatile Multi Spring Adjustability (VMSA) and visual feedback through the LED lighted control panel



#### **Comfortable harness**

The comfortable harness is especially made for Lode to guarantee a comfortable fitting for the patient. The patient will not be burdened by the compensation weight, enabling the most comfortable way of rehabilitation. The harness is available in various sizes: small, standard and extra large.

#### **Single Point Suspension**

to move

Thanks to the single point suspension dynamic up and down movement is possible and allows the patient to move and rotate freely.



#### Body Weight Supported Treadmill Training

Patients that have to learn how to walk again can be anxious. The Lode BWSS offers various tools to take away this anxiety. Of course, safety is one of the most important aspects that Lode takes into consideration. That is why a fallstop comes standard with this system. In case a patient impends to fall, the fallstop will make sure that the patient will remain in the system and will not be injured.

Many patients that have to learn how to walk again are mobile through a wheelchair. The Lode BWSS takes exactly these patients into full consideration: with the optional entrance plate the patient can be placed on the treadmill in the wheelchair. The harness can easily be put on to the sitting patient. Then the patient can stand up directly from the wheelchair with support of the Lode BWSS thanks to the lifing function. This function is suitable for patients up to 160 kg. Now the patient can start the therapy without much effort.






### XL

#### Reliable support system for rehabilitation

Body Weight Support System XL can a.o be extended with the following options:





### XL

Accuracy

#### Reliable support system for rehabilitation

### Specifications



Weight support steps	0.1 kg	0.2 lbs
Comfort		
Versatile Multi-Spring Adjustability (VSMA) Min.	0 kg	
Versatile Multi-Spring Adjustability (VSMA) Max.	80 kg	
Allowed user weight	160 kg	352.7 lbs
Maximum lifting weight	160 kg	352.7 lbs
Dimensions		
Product length (cm)	160 cm	63 inch
Product width (cm)	225 cm	88.6 inch
Product height	250 cm	98.4 inch
Product weight	110 kg	242.5 lbs
Inner width	193 cm	76 inch
Inner width with optional handrails	149 cm	58.7 inch
Maximum patient height	215 cm	84.6 inch

#### Order info

Partnumber:

932832

\*Specifications are subject to change without notice.



Distributed by Haganjordet 30 1351 Rud Norway

Tel:

Lode B.V. Zernikepark 16 9747 AN Groningen The Netherlands Tel: +31 50 5712811 Fax: +31 50 5716746 E-mail: ask@lode.nl Internet: www.lode.nl